



# THULI TABLES, INC.

800.458.4854 608.935.9300 fax 608.935.9192 office@thulitables.com www.thuli.com

## OWNER'S MANUAL ADDENDUM

## 300 Abdominal Swing-Away



This addendum to the 300 Owner's Manual will familiarize you with the features and operation of the Abdominal Swing-Away and Interchangeable Chest Cushions.

## Abdominal Swing-Away Section: Operation

### ABDOMINAL SWING-AWAY: POSITIONING



Activate the lever on either side of the table while depressing the abdominal section. This may be done before the patient lays on the table or while the patient is on the table. Adjust the position of the abdominal swing-away to the patient's comfort.



This abdominal swing-away section offers prone comfort to pregnant patients or to those of us with “a little extra around the middle.”

## Interchangeable Chest Cushions

### STANDARD CHEST CUSHION:

This cushion has a foam density that matches the firmness of the rest of the table. This is likely the cushion you will use for most patients.

### “FEMALE FRIENDLY” CHEST CUSHION:

This cushion provides a combination of soft foam in the breast area and firmer foam surrounding, insuring adequate thoracic support when the patient is lying prone. The firmer sternal cushion allows for performing “anterior thoracic” adjustments.

### CHEST CUSHION: PLACEMENT



The chest cushions are attached by way of four magnets. To remove the cushion, simply grasp at both sides and twist toward the headpiece with a lifting motion. The cushion is held in correct position by alignment of two holes on the underside of the chest cushion with the two nickel plated acorn nuts on the thoracic drop plate.

### CHEST CUSHION: STORAGE



When not using the chest cushion, consider:

1. Setting it aside.
2. Storing it in the space vacated by the extended ankle rest section (photo above).
3. Use it as a bolster under the patient's ankles when lying prone.
4. Use it as a bolster under the patient's knees when lying supine.