

# THULI ELEVATION TABLES

## MOTOR CHECK LIST


---

### MOTOR DOES NOT WORK


1. The outlet may not be working or the breaker is tripped. Try plugging in another device.
2. The electrical connection in the back of the motor may be disengaged.
3. The air lines from the foot switch to the motor may be not be fully connected or they may be kinked or broken. If they are kinked, try to straighten them. If broken, we sell replacement air lines.
4. The foot switch may be worn out. Uplug the air lines from the foot switch and blow into both of them, one at a time. If the foot switch is worn out, the motor will operate by doing this. We sell replacement foot pedals.
5. The wires inside the plug may be loose or disconnected.  
To check this, first unplug the motor from the electricity. Follow the steps below.



Remove the plug from the motor by lifting up on the tab (using a flat screwdriver).



Remove the cover by removing the 2 screws (using a Philips screwdriver).



Loosen the 3 screws (using a flat screwdriver) and re-connect the wires as follows:

**110 Motors:**

White wire goes into L1 (left)  
Green wire goes into middle  
Black wire goes into N (right)

**220 Motors:**

Brown wire goes into L1 (left)  
Green/Yellow goes into middle  
Blue goes into N (right)

6. If none of the above, the motor must be sent in for repair. Please contact us.

### MOTOR IS MAKING NOISE

1. Squealing or screeching noise:  
The electronic brake has likely failed and must be sent in for repair.
2. Clicking noise:  
This is normal and may get louder with age, but functionally the motor is fine.

### MOTOR / TABLE MOVEMENT

When the table is in an elevated position, some movement of the table is normal. If the movement becomes excessive and interferes with function, it is likely that the motor needs replacing. Please contact us.



800.458.4854  
608.935.9300  
www.thulitables.com  
office.thulitables.com