


**QUICK GUIDE: OPENING UP THE SPORT TABLE**

1. Lay the table onto its contoured thoracic section.
2. If your table has an ankle rest extension (optional), loosen the locking knob and pull out the ankle rest.
3. Release the eye hook from the black metal bracket located between the headpiece cushions.

**IMPORTANT NOTE: RAISING THE FRONT LEGS**

Swing the headpiece and front set of legs up simultaneously by grasping the front legs with one hand and the headpiece with the other hand. See owner's manual pg. 6.



**QUICK GUIDE: OPENING UP THE SPORT TABLE**

1. Lay the table onto its contoured thoracic section.
2. If your table has an ankle rest extension (optional), loosen the locking knob and pull out the ankle rest.
3. Release the eye hook from the black metal bracket located between the headpiece cushions.

**IMPORTANT NOTE: RAISING THE FRONT LEGS**

Swing the headpiece and front set of legs up simultaneously by grasping the front legs with one hand and the headpiece with the other hand. See owner's manual pg. 6.