

Side Arm Rests



Do you use sometimes need a wider table for larger patients, Activator, Craniosacral or Soft Tissue therapies? Make it easier on yourself and more comfortable for the patient with added width in the upper torso area. Attach quickly. Available only on specific Thuli Tables (see below). Set of 2, each 4" x 26".

Tour Portables, 300 Stationary & 300 Oak Elevation Tables
Increases 20" table width from to 28" - 32" (adjustable)

Tour on Conversion Base & 300 Swing-Away Table
Increases 20" table width from to 28" (not adjustable)

PLEASE NOTE:

- The Side Arm Rests must be ordered at the time of table for Thuli Table installation.
- The Arm Rest on the right side must be removed when using the Thoracic and Pelvic Drops.
- They must be carried separately from the Tour Portable during transport.



THULI TABLES

800.458.4854 / 608.935.9300

608.935.9192 Fax

www.thulitables.com

office@thulitables.com

facebook.com/thulitables

youtube.com/thulitables