

Sport Flexion Distraction



Add Flexion Distraction to your Sport Portable with this ingenious, mechanical device. Simply fold the table's back legs into their storage position and effortlessly attach the Device in seconds. It also quickly disassembles for portability so you can take it with you. Component parts are individually protected in a custom Travel Case (optional).

Lateral Flexion: Available is only on the Tour Portable.

RETROFIT

Available for 18" - 24" Sports

Retrofit Kit needed for Sports built before July 2013

FLEXION DISTRACTION OPTIONS

Star Base: Black (Red Oak or Walnut optional)

Operating Knob: Black, Red Oak or Walnut

ACCESSORIES

Flexion Distraction Travel Case

Comfort Pillow for proper patient positioning

Plastic floor mat for optimum rolling

FLEXION DISTRACTION DEVICE

Flexion Range: 14 - 17" Sport height: Not available
18" Sport height: 0 - 10 degrees
19" - 21" Sport height: 0 - 11.5 degrees
22" - 24" Sport height: 0 - 14.5 degrees
Locks in position at any angle

Star Base: Black (Red Oak or Walnut optional), 21" diameter

Operating Mechanism: Gas spring (smooth & quiet operation)

Operating Knob: Red Oak, Walnut or Black. Ergonomic hand-turned hardwood knob, 3.25" diameter

Wheels: 5 specialty wheels

Leg Weights: (2) 10 lb counterweights provided

Conversion Time: Attaches to or detaches from the Sport in one minute

Retrofit Kit: Available for 18" - 24" Sports prior to July 2013

Weight: 37 lbs

Ship Weight: 42 lbs

62 lbs (with leg weights)

68 lbs (with leg weights & Travel Case)

Box size: 24" x 24" x 8" (without Travel Case)

24" x 24" x 10" (with Travel Case)



THULI TABLES

800.454.4858 / 608.935.9300

608.935.9192 Fax

www.thulitables.com

office@thulitables.com

facebook.com/thulitables

youtube.com/thulitables