

# Tour Flexion Distraction and Lateral Flexion

---



## FLEXION DISTRACTION DEVICE

Add Flexion Distraction to your Tour Portable with this ingenious, mechanical device. Simply fold the table's back legs into their storage position and effortlessly attach the Device in seconds. It also quickly disassembles for portability so you can take it on the road. Component parts are individually protected in a custom Travel Case (optional).

## LATERAL FLEXION BASE

Must be used in conjunction with the Flexion Distraction Device. The Lateral Flexion Base offers 30 degrees right and 30 degrees left flexion and can be locked in at any position along the 60 degree arc. Multiple, compound angles are also available. Speciality wheels offer quiet, smooth and instant directional change. We do not consider the Lateral Flexion Base to be "portable".

## RETROFIT

Available for 20"- 24" Tours built after 2001  
Retrofit Kit needed for tables 2001 - March 2013

## FLEXION DISTRACTION OPTIONS

Star Base: Black (Red Oak or Walnut upgrade)  
Operating Knob: Black, Red Oak or Walnut

## ACCESSORIES

Flexion Distraction Travel Case  
Comfort Pillow for proper patient positioning  
Plastic floor mat for optimum rolling

## FLEXION DISTRACTION DEVICE

---

**Flexion Range:** 17" Tour height: Not available  
18" Tour height: 0 - 11 degrees  
19"- 21" Tour height: 0 - 12 degrees  
22"- 24" Tour height: 0 - 15 degrees  
Locks in position at any angle

**Star Base:** Black (Red Oak or Walnut optional), 21" diameter

**Operating Mechanism:** Gas spring (smooth & quiet operation)

**Operating Knob:** Red Oak, Walnut or Black. Ergonomic hand-turned hardwood knob, 3.25" diameter

**Wheels:** 5 wheels (smooth & quiet operation)

**Leg Weights:** (2) 10 lb counterweights provided

**Conversion Time:** Attaches to or detaches from Tour in one minute

**Retrofit Kit:** Available for 18" - 24" Tours from 2001 - May 2013

**Device Weight:** 37 lbs

**Box size:** 24" x 24" x 8" (without Travel Case)

24" x 24" x 10" (with Travel Case)

**Shipping Weight:** 38 lbs

58 lbs (with leg weights)

64 lbs (with leg weights & Travel Case)

## LATERAL FLEXION BASE

---

**Lateral Range:** 17" - 19" Tour height: Not available  
20" - 24" Tour height: 60 degree lateral range  
(30 degrees right, 30 degrees left)  
Locks in position at any angle

**Base:** 40" x 20" x (4" less than table height)

Red Oak shroud over steel base, non-skid feet

**Conversion Time:** Attaches to and detaches from Tour in minutes

**Weight:** 57 lbs

**Ship Weight:** 83 lbs

**Box size:** 57" x 26" x 7"